

OHAUS Welcome To Ohaus




การเลือกใช้เครื่องชั่งให้เหมาะกับงาน
How to choose the right balance

OHAUS SI units? Basic Training

■ Mass (kg)

Structure of the SI units.



length **m**


Mass **kg**

Time **s**

Current **A**

Temp **K**


6/5/2013

 **Basic Training**

What is Weight? Mass vs. Weight


Mass = Constant unit of the amount of matter an object possesses. Stays the same no matter where the measurement is taken.

Mass = Volume x Density



Typical mass measurement unit is kg or gram

6/5/2013 3

 **Basic Training**


What is Weight? Mass vs Weight

Weight = Heaviness of an object. Depends on the gravity exerted on the object multiplied by its mass.

Weight = Mass x Gravity

Actual weight measurement unit is the Newton

6/5/2013 2




What is Weight?


Basic Training

“A force whose dimension depends on gravity”

Unit: 1N (Newton)



The weight of a stone on the Moon is lower than the weight of the same stone on Earth. However, its mass remains the same.




The display on scales indicates the weight of the object weighed in its current field of gravity, expressed in grams, i.e. the “apparent mass”.

Given that weight acceleration g varies between locations, balances must be adjusted if moved.

When transported, the weight displayed by the balance changes by approx. 62 mg per kg and degree of latitude

5

6/5/2013



How Much is 1 Gram?


Basic Training

Digits	Readability (g)	Readability (mg/ µg)	Type of balance
0	1 g	1000 mg	Precision balance
1	0.1 g	100 mg	Precision balance
2	0.01 g	10 mg	Precision balance
3	0.001 g	1 mg	Precision balance (but with a draft shield)
4	0.0001 g	0.1 mg	Analytical balance
5	0.00001 g	0.01 mg	Semi-Micro balance
6	0.000001 g	1 µg	Micro balance
7	0.0000001 g	0.1 µg	Ultra-Micro balance

1 t	=	1000 kg
1 kg	=	1000 g
1 g	=	1000 mg
1 mg	=	1000 µg

6

6/5/2013



Basic Training


Scales or Balances

Both are devices that measure weight by mechanical or electronic means.

- In general a balance is a device with a readability of 0.1g or better. Scales have readabilities of 1g or higher.
- Determined by convention usually lab balance as opposed to lab scale or truck scale as opposed to truck balance.

7

6/5/2013



Basic Training

Capacity and Readability

Capacity


- The total weight that can be measured on a balance.

Readability

- The smallest weight that can be displayed on a balance.

8

6/5/2013



Basic Training

Range and Units


Weighing Range

- The minimum load value to the capacity of a balance.

Weighing Units

- The units of measure a balance will display.
 - Grams or Kilograms are usually the default, other units may have be turned on before they can be used.

6/5/2013 9




Basic Training

Resolution

Determined by dividing the capacity of a balance by its readability

- A balance with 1000g capacity and .01g readability would have a resolution of 10,000 counts
- Also referred to as scale divisions (d).
- $d = 0.01g$

6/5/2013 10



Basic Training

OHAUS

Criteria for purchasing Lab Balance

- Capacity** – What is the most it will weigh?
- Readability** – What is the least it will weigh?
- Linearity** – What features does it need?
- Repeatability** – Where will it be used?
- Minimum Weight** – How reliable is it?
- Price** – How expensive is it?

12

6/5/2013

OHAUS What is the right balance? Basic Training


We would like to produce 2 kg emulsion in a beaker weighing 500g. We need the following ingredients:

Water :	70%
Oil:	20%
Emulgator:	5%
Moisturizer:	3%
Stabilisator:	1.8%
Fragrance:	0.2%

→ 4 g

Each component shall be weighed with an accuracy of at least 1%.
If possible, we would like to work on a single balance.

Which balance do you recommend ?



6/5/2013 13

OHAUS What is the right balance? Basic Training




Can I use a balance with a readability of 1 g to weigh my 4g of fragrance with an accuracy of 1%?

So I need a balance with a higher readability, such as 0.1g, or 0.01g or maybe even 0.001g?

The readability of a balance can be set to any value electronically – which does NOT reflect the accuracy of the instrument.
What we need to know is the **measurement uncertainty** – and as its most important consequence - the **minimum weight** of a balance.

6/5/2013 14




Basic Training

Accuracy

The accuracy of the measurement refers to how close the measured value is to the true or accepted value.

15

6/5/2013




Basic Training

Factors affecting accuracy

- Temperature
- Air currents
- Calibration accuracy
- Vibration
- Changes In temperature
- Changes in atmospheric pressure
- Method of weighing

16

6/5/2013




Basic Training

Accuracy

Can only be measured when the balance is set up where it is going to be used using the method that is going to be used.

17

6/5/2013



Basic Training


Repeatability

How close are successive readings of the same unknown using the same method.

- Expressed as Standard Deviation

18

6/5/2013



Basic Training

Linearity


Linearity of a balance is the measure of the error between the displayed weight reading and the true value of the weight placed on the balance at different points throughout the range of the balance.

19

6/5/2013


What is the right balance
Basic Training

EX4202
Readability = 0.01g




$$U [g] = 0.012 + 3.5e-6 \times \text{Weight}$$

Load on balance [g]	Absolute Measurement Uncertainty [mg]
0.01	12
0.1	12
1	12
10	12
100	12
1000	17
4100	26

20

6/5/2013




What is the right balance

Basic Training

$U [g] = 0.012 + 3.5e-6 \times \text{Weight}$


EX 4202
Readability = 0.01g



Load on balance [g]	Absolute Measurement Uncertainty [mg]	Relative Measurement Uncertainty [%]
0.01	12	
0.1	12	
1	12	
10	12	
100	12	
1000	17	
4100	26	0.0006

21 $\frac{26\text{mg}}{4100\text{g}} = 0.0006\%$

6/5/2013




What is the right balance

Basic Training


$U [g] = 0.012 + 3.5e-6 \times \text{Weight}$

EX 4202
Readability = 0.01g





Load on balance [g]	Absolute Measurement Uncertainty [mg]	Relative Measurement Uncertainty [%]
0.01	12	120
0.1	12	12
1	12	1.2
10	12	0.12
100	12	0.012
1000	17	0.0017
4100	26	0.0006

±12 %
88 mg ≤ Weight ≤ 112 mg



6/5/2013


What is the right balance
Basic Training



Load on balance [g]	Absolute Measurement Uncertainty [mg]	Relative Measurement Uncertainty [%]
0.01	12	120
0.1	12	12
1	12	1.2
10	12	0.12
100	12	0.012
1000	17	0.0017
4100	26	0.0006

not OK

← Limit ?

OK

6/5/2013
Z3



Basic Training

การเลือกซื้อเครื่องชั่งอุตสาหกรรม





6/5/2013


Basic Training

ทำไมต้องเลือก

1. ต้องการให้ได้เครื่องชั่งที่เหมาะสมกับการใช้งานมากที่สุด
 - ความถูกต้องในการชั่ง (Accuracy), ลักษณะการใช้งาน และ อื่นๆ
2. เพิ่มประสิทธิภาพในขบวนการผลิต
3. เพื่อรักษาให้เครื่องชั่งอยู่กับเราได้นานที่สุด

ถ้าไม่ถูกต้อง

1. ขบวนการผลิตขาดประสิทธิภาพ
2. องค์กรเสียหาย
3. เสียง่าย เกิดปัญหาฉุกเฉิน

6/5/2013


Basic Training

สิ่งที่เราต้องทราบ


1. พิกัดเครื่องชั่ง
2. ค่าอ่านละเอียด
3. ขนาดแทนชั่ง
4. โครงสร้างของเครื่องชั่ง
5. ชนิดของโหลดเซลล์
6. การประยุกต์ใช้งาน








6/5/2013


Basic Training

1. พิกัดเครื่องชั่ง Weighing Capacity

- พิกัดสูงสุดของเครื่องชั่งหาจาก
 - น้ำหนักสูงสุดที่ต้องการชั่ง ÷ 0.7
 - เครื่องชั่งมีความสามารถสูงสุดที่บริเวณ 70% ของพิกัด



ตัวอย่าง ต้องการชั่งสินค้า 30 กก.

พิกัดสูงสุด = $30/0.7 = 42$ กก.

เลือกพิกัด = 60 กก.

6/5/2013


Basic Training

สิ่งที่เราต้องทราบ

1. พิกัดเครื่องชั่ง
2. ค่าอ่านละเอียด
3. ขนาดแท่นชั่ง
4. โครงสร้างของเครื่องชั่ง
5. ชนิดของโพลดเซลล์
6. การประยุกต์ใช้งาน





6/5/2013


 **Basic Training**


2. คำอ่านละเอียด

Readability

- คำอ่านละเอียดหาได้จาก
 - ความผิดพลาดสูงสุดของสินค้าที่ต้องการตรวจสอบ
 - ISO10012 : Quality Assurance Requirement for Measuring Equipment-
Part 1: Metrological Confirmation System for Measuring Equipment
4.3 Confirmation System
GUIDANCE
"The error attributable to calibration should be as small as possible. In most areas of measurement, it should be no more than one third and preferably one tenth of the permissible error of the confirmed equipment when in use".

6/5/2013

 **Basic Training**




ผิดพลาด = 1 ใน 3

ผิดพลาด = 1 ใน 3

ผิดพลาด = 1 ใน 3

6/5/2013


Basic Training

2.ค่าอ่านละเอียด (Readability)

- ค่าอ่านละเอียดหาได้จาก
 - เครื่องชั่งที่ใช้ตรวจสอบสินค้าควรมีความผิดพลาด $\frac{1}{3}$ หรือ $\frac{1}{10}$ ของความผิดพลาดของสินค้า
ตัวอย่าง ต้องการชั่งสินค้า 30 กก

#1 ความผิดพลาดไม่เกิน	+/- 1%	= +/- 300 g
เครื่องชั่งควรมีความผิดพลาดไม่เกิน	= 300/3	= +/- 100 g
#2 ความผิดพลาดไม่เกิน	+/- 0.1%	= +/- 30 g
เครื่องชั่งควรมีความผิดพลาดไม่เกิน	= 30/3	= +/- 10 g

6/5/2013


Basic Training




2.ค่าอ่านละเอียด Readability

- ค่าอ่านละเอียดหาได้จาก
 - นำค่าความผิดพลาดเครื่องชั่งหารด้วย 3 (Ref. to MPE chart)
ตัวอย่าง ต้องการชั่งสินค้า 30 กก

#1 ความผิดพลาดไม่เกิน	+/- 1%	= +/- 300 g
เครื่องชั่งควรมีความผิดพลาดไม่เกิน	= 300/3	= +/- 100 g
ค่าอ่านละเอียดเครื่องชั่ง	= 100g/3 = 33 g	= 20g
#2 ความผิดพลาดไม่เกิน	+/- 0.1%	= +/- 30 g
เครื่องชั่งควรมีความผิดพลาดไม่เกิน	= 30/3	= +/- 10 g
ค่าอ่านละเอียดเครื่องชั่ง	= 10g/3 = 3 g	= 2 g

6/5/2013

OHAUS Basic Training

สิ่งที่เราต้องทราบ

1. พิกัดเครื่องชั่ง
2. ค่าอ่านละเอียด
3. ขนาดแท่นชั่ง
4. โครงสร้างของเครื่องชั่ง
5. ชนิดของโหลดเซลล์
6. การประยุกต์ใช้งาน



6/5/2013

OHAUS Basic Training

3. ขนาดแท่นชั่ง Platform Size

- สิ่งของที่ต้องการชั่งจะต้องอยู่ในพื้นที่ของแท่นชั่ง



6/5/2013


Basic Training

สิ่งที่เราต้องทราบ

1. พิกัดเครื่องชั่ง
2. คำอ่านละเอียด
3. ขนาดแท่นชั่ง
4. โครงสร้างของเครื่องชั่ง
5. ชนิดของโหลดเซลล์
6. การประยุกต์ใช้งาน

6/5/2013


Basic Training

4. โครงสร้างของเครื่องชั่ง

Platform & Terminal Structure (SS304,316,MS)

- พิจารณาจากวัสดุที่ต้องการชั่งน้ำหนัก (กรด, ด่าง)
- พิจารณาจากสภาพแวดล้อมการใช้งาน (Hazardous Area, wash down, Dust, etc.)



6/5/2013

OHAUS Basic Training

4. โครงสร้างของเครื่องชั่ง
Platform & Terminal Structure (SS304,316,MS)



MS Structure



SS Structure

6/5/2013

OHAUS Basic Training

การป้องกันน้ำและฝุ่นละออง IP Protection IP "XX"
IEC Standard

Dust Protection


IP Xx
IP 6x

Water Protection

IP xX
IP x5

~~IP 65~~


6/5/2013


Basic Training

Degrees of Protection (IP Standard)

การป้องกันของแข็ง (Degrees of Protection) และ IP ตามมาตรฐาน DIN40050/1980 และ IEC 529		การป้องกันของแข็ง (Degrees of Protection) และ IP ตามมาตรฐาน DIN40050/1980 และ IEC 529		การป้องกันของแข็ง (Degrees of Protection) และ IP ตามมาตรฐาน DIN40050/1980 และ IEC 529	
IP	รายละเอียด	IP	รายละเอียด	IP	รายละเอียด
0	ไม่มีการป้องกัน	0	ไม่มีการป้องกัน	0	ไม่มีการป้องกัน
1	สามารถป้องกันของแข็งที่ มีเส้นผ่านศูนย์กลางมากกว่า หรือเท่ากับ 50 มม. ที่มาจากทุกทิศทาง	1	สามารถป้องกันน้ำที่ตกลงมาในแนวตั้งได้	1	สามารถป้องกันแรงกดกรรพรม ของวัตถุที่มีน้ำหนัก 150 กรัม ที่ปล่อยมาจากที่สูง 1.5 ซม. (การกรรพรมของพลังงาน 0.25 จูล)
2	สามารถป้องกันของแข็งที่ มีเส้นผ่านศูนย์กลางมากกว่า หรือเท่ากับ 12 มม. ที่มาจากทุกทิศทาง	2	สามารถป้องกันน้ำที่ตกลงมาใน แนวตั้ง แต่ในแนวเอียง 15 องศา กับแนวตั้ง	2	สามารถป้องกันแรงกดกรรพรม ของวัตถุที่มีน้ำหนัก 250 กรัม ที่ปล่อยมาจากที่สูง 1 ซม. (การกรรพรมของพลังงาน 0.25 จูล)
3	สามารถป้องกันของแข็งที่ มีเส้นผ่านศูนย์กลางมากกว่า หรือเท่ากับ 2.5 มม. ที่มาจากทุกทิศทาง	3	สามารถป้องกันน้ำฝนที่ตกลงมา ได้โดยพายุที่อาจตกลงมาใน แนวตั้งมุม 60 องศา กับแนวตั้ง	3	สามารถป้องกันแรงกดกรรพรม ของวัตถุที่มีน้ำหนัก 500 กรัม ที่ปล่อยมาจากที่สูง 20 ซม. (การกรรพรมของพลังงาน 0.25 จูล)
4	สามารถป้องกันของแข็งที่ มีเส้นผ่านศูนย์กลางมากกว่า 1 มม. ที่มาจากทุกทิศทาง	4	สามารถป้องกันหยดน้ำ หรือ น้ำสาด ที่มาจากทุกทิศทาง	4	สามารถป้องกันแรงกดกรรพรม ของวัตถุที่มีน้ำหนัก 500 กรัม ที่ปล่อยมาจากที่สูง 40 ซม. (การกรรพรมของพลังงาน 0.25 จูล)
5	สามารถป้องกันฝุ่นได้	5	สามารถป้องกันน้ำที่ฉีดออกมาจาก หัวฉีดในทิศทางต่างๆ	5	สามารถป้องกันแรงกดกรรพรม ของวัตถุที่มีน้ำหนัก 1.5 กก. ที่ปล่อยมาจากที่สูง 40 ซม. (การกรรพรมของพลังงาน 0.25 จูล)
6	สามารถป้องกันฝุ่นได้อย่างสมบูรณ์	6	สามารถป้องกันคลื่นตามยาวที่เกิดจาก คลื่นวิทยุแรง และ การฉีดน้ำอย่างแรง	6	สามารถป้องกันแรงกดกรรพรม ของวัตถุที่มีน้ำหนัก 5 กก. ที่ปล่อยมาจากที่สูง 40 ซม. (การกรรพรมของพลังงาน 0.25 จูล)
		7	สามารถป้องกันอันตรายที่เกิดจาก น้ำท่วมชั่วคราว		
		8	สามารถป้องกันอันตรายที่เกิดจาก น้ำท่วมถาวรได้		

6/5/2013


Basic Training

สิ่งที่เราต้องทราบ

1. พิกัดเครื่องขั้ว
2. ค่าอ่านละเอียด
3. ขนาดแท่นขั้ว
4. โครงสร้างของเครื่องขั้ว
5. ชนิดของ โหลดเซลล์
6. การประยุกต์ใช้งาน

6/5/2013


Basic Training

5.ชนิดของ โหลดเซล Load Cell

- พิจารณาจากวัสดุที่ต้องการชั่งน้ำหนัก (กรด, ด่าง)
- พิจารณาจากสภาพแวดล้อมการใช้งาน (Hazardous Area, wash down , Dust ,etc.)
- พิจารณาจากลักษณะการใช้งาน (Heavy Duty)

6/5/2013


Basic Training

5.ชนิดของ โหลดเซล Load Cell

วัสดุโหลดเซล

- อลูมิเนียม Aluminum Alloy
 - ไม่ทนต่อการกระแทก Impact
- เหล็กกล้าผสมต่ำ Low Alloy Steel
 - แข็งแรง ทนทานต่อการกระแทกได้ดีกว่า
 - ไม่ทนต่อการกัดกร่อน Corrosion Resistance
- เหล็กกล้าสแตนเลส Stainless Steel
 - แข็งแรง ทนทาน
 - ทนการกัดกร่อน Corrosion Resistance



6/5/2013


Basic Training

สิ่งที่เราต้องทราบ

1. พิกัดเครื่องชั่ง
2. คำอ่านละเอียด
3. ขนาดแท่นชั่ง
4. โครงสร้างของเครื่องชั่ง
5. ชนิดของโหลดเซลล์
6. การประยุกต์ใช้งาน





6/5/2013


Basic Training

Portable Scale

- ลักษณะการใช้งาน
 - ตั้งบนโต๊ะ
- พิกัดน้ำหนัก(พิกัดน้ำหนักไม่มาก)
 - 3 kg.-30kg, Platform 19x25 cm.

Bench Scale

- ลักษณะการใช้งาน
 - วางกับพื้น
 - หัวกับแท่นแยกกันเพื่อความสะดวก หรือป้องกันการเปื้อน
- พิกัดน้ำหนัก
 - 15 kg-300 kg, Platform 30x30 cm – 61 x61 cm.





6/5/2013


Basic Training

Floor Scale

- ลักษณะการใช้งาน
 - วางกับพื้น . ผึ่งบ่อ
 - หัวกับแท่นแยกกัน
 - สำหรับงานที่มีการกระแทกและของมีขนาดใหญ่
- พิกัดน้ำหนัก
 - 1500 kg.-3,000 kg,
 - Platform 1mX1m-1.5mX1.5m





6/5/2013



Thank You