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- Source waters (environmental testing)
 - Drinking water
 - Ground Water
 - Surface Water
- Reversed Osmosis Water (RO-water)/Cooling water
- Efficiency of water polishing plant (Control)
- High Purity Water or Purified Water (PW)
- Pharmaceutical Grade Water, WFI (Water For Injection)
- Clean In Place (CIP)

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
Oxidation Techniques

- Combustion
- Heated Persulfate
- Ultraviolet & Persulfate
- Ultraviolet and Titanium dioxide
- Ultraviolet

Basic Building Blocks of TOC Analysis

1. Syringe
2. Process / On-line
3. Sample loop
4. Autosampler
5. Solids option

1. Combustion
2. UV/Persulfate
3. Persulfate at 100°C
4. UV Only
5. UV / TiO₂

1. Non-dispersive Infrared Detector
2. Conductivity
3. Coulometric

1. Display
2. Printout
3. Computer

Combustion Oxidation

High temperature furnace + air/O₂ + Pt catalyst

Advantages

High oxidation efficiency



Oxidize some particles

Disadvantages

Catalyst lifetime

Reagents required

Heated Persulfate Oxidation

Advantages

High oxidation efficiency

No catalyst required

Appropriate – Cleaning Validation

Disadvantages

Reagents required

Ultraviolet

Advantages

No reagents

Disadvantages

Low oxidation efficiency

Can not oxidize particles

Lamp replacement

Subtraction for TOC

UV Persulfate Oxidation

Advantages

Expensive Reagents not required

Disadvantages

Low Oxidation Efficiency

Required pre-screening of samples

Lamp replacement – Expensive

Does not oxidize particles

Oxidation Techniques

<i>Technique</i>	<i>Range</i>	<i>Oxidation Efficiency</i>
Combustion	500 ppb – 5000 ppm	very high
Heated Persulfate	2 ppb – 125 ppm (automatic)	high
	2 ppb – 10% (syringe)	high
Ultraviolet & Persulfate	<1 ppb – 1 ppm (unscreened)	low
	< 1 ppb – 40 ppm (prescreened)	low
Ultraviolet /TiO ₂	2 ppb – 1 ppm	low
Ultraviolet	5 ppb – 1 ppm	very low

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

- Good oxidation of organics
- Low and stable system blanks
- Relatively fast sample analysis time
- Dissolved organic carbon analysis
- Particle free (filtered) (or suspended matter <0.15 mm ID (USEPA 415.1))

Heated Persulfate Oxidation

- Easy to validate (QA/QC control (21 CFR part 11))
- High accuracy & precision
- Measuring in inert gas (no interferences) by NDIR
- No cross-contamination
- Low cost of operation
- Lower detection limits
- Lower maintenance
- TOC ppb levels to low ppm, but can be up to 30,000 ppm

Here precision and accuracy at low ppb levels is important

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$$\text{TC} = \text{TIC} + \text{TOC}$$

Combustion Oxidation

- Good Oxidation of Organics
- Suitable for Dirty Water Samples
- Good Precision & Accuracy at ppm levels
- Good at oxidizing these material, especially hard to oxidize compounds such as leaf material, polystyrene, etc.

Combustion versus Persulfate

Combustion	Persulfate
Dirty Water Samples	Clean Water Samples
High detection limits (ppm)	Lower detection limits (ppb)
Good for difficult to oxidize compound or high molecular weight organics	Used in analysis of drinking water and pharmaceutical water

Lab Technique Vs. On-line Technique

Lab

- Easy to validate (QA/QC control)
- delayed data
- TOC as NPOC
- High Accuracy & Precision (large sample size (0.5-20 ml))
- Measuring in inert gas (No interferences) by NDIR
- Priority sampling
- Fully self serviceable!
- Wide analytical range
- Source, RO, PW, WFI, CIP
- No cross-contamination

On-line

- More difficult to validate (process environment and autom. (PC or ...))
- Real time data
- TOC as TC-IC
- Good precision but worse accuracy (low sample size (nl))
- Interference from Ionic species & Halogenated Org. (Na, K, Cl, O3, NH3...) by Conductivity based detection
- No priority sampling
- Narrow analytical range
- Sensitive for cross-contamination

Reagents used in the analysis

Reagents for the Aurora consist of

- reagent water (recommends using distilled or deionized water containing less than 200 ppb TOC)
- phosphoric acid required (5% Phosphoric acid solution)
- sodium persulfate ($\text{Na}_2\text{S}_2\text{O}_8$) (required for analysis is a 10% solution, or 100 g/L sodium persulfate)

Examples of Applications

Typical Levels of Carbon in Water

Type of Water	TOC (mg/L)	IC (mg/L)
Groundwater	<1	5–300
Surface Water	<10	~200
Surface Water, Polluted	>10	~200
Process Water (Ultra-High Purity)	0.005–0.1	<0.05
Drinking Water	1–5	~50
Water for Injection (WFI)	0.05–0.5	<1
Clean in Place (CIP)	0.20–50 ppm	<1


Ultrapure (Distilled) Water

NPOC only Internal					
<i>TIC Sparge</i>	2:00		<i>Sample Vol.</i>	8mL	
<i>TOC React</i>	2:00		<i>Acid Vol.</i>	0.5mL	
<i>TOC Detect</i>	3:00		<i>Persulfate Vol.</i>	1mL	
20psi EPC standby and detect					
120mL/min EFC					
RW subtraction					
RW	866	5.40		RF	0.3639
1ppm	22921	0.37		R2	1.000
5ppm	110581	0.32		Offset	887
10ppm	220905	0.10		RB	752

Ultrapure (Distilled) Water

<u>Distilled Water</u>			
<i>Results</i>			
	ppm	%RSD	ppb error
1	0.049	8.10	3.97
2	0.045	2.56	1.15
3	0.043	2.71	1.17
4	0.042	4.46	1.87
5	0.043	3.09	1.33
6	0.040	1.83	0.73
7	0.040	6.67	2.67
8	0.044	1.58	0.70
9	0.041	4.55	1.87
10	0.042	2.02	0.85
11	0.041	2.53	1.04
12	0.039	3.22	1.26
13	0.041	1.74	0.71
14	0.042	3.48	1.46
15	0.041	3.14	1.29
16	0.043	3.36	1.44
17	0.042	2.59	1.09
18	0.042	2.86	1.20
19	0.042	1.34	0.56
20	0.042	2.93	1.23
Average (omitting #1)	0.042	2.98	1.24
incl #1	0.042	3.24	1.38

Drinking Water: Continuous Monitoring

Thank You